



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Capillary Tube Clamp

## Product Image

---



## Description

---

The apparatus comprises a metal frame, on which four capillary tubes of 10cm long and of different internal diameters are clamped with the help of small metal brackets. The 9mm rod is also attached to metal plate, so that complete clamp can be held in any stand for determining the surface tension of liquid by capillary rise method.

Capillary Tube Clamp

**Catalog No. 100406**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.