



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cell Wash Centrifuge

Product Image



Description

Model : T-5

Specifications:

- Full body made of ABS plastic.
- Fitted with extra long life motor without carbon.
- Programmable speed regulator to control speed in range of 500, 1000, 1500, 2000, 2500, 2850 R.P.M.
- 0 to 99 minutes timer which can be set in interval of 1 minute.
- Supplied with Swing out Rotor to hold 2 plates of 96 wells each.

Physical Parameters:

Net weight	12.5 Kg. Approx.
Gross weight with box	14 Kg.
Dimensions of machine	L = 395mm, W = 360 mm, H = 280 mm
Dimensions of outer box	510 x 510 x 390 mm

-
Cell Wash Centrifuge

Catalog No. 100761

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.