



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Centrifuge Manual (Hand Type)

## Product Image

---



## Description

---

For 15ml. test tubes with swing out head and polyamide buckets to house test tubes. It has clamping device for fixing on table. It gives 1500 to 2500 revolutions per minute.

<b>Catalog No.</b>	<b>Particulars</b>
JZ00540	Centrifuge Manual (Hand Type), 2 tubes
JZ00541	Centrifuge Manual (Hand Type), 4 tubes

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.