



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Centrifuge Tube, Micro, Plastic (Dnase/Rnase Free Certified)

## Product Image



## Description

---

**Material:** Polypropylene

Theae tubes assures good transparency and clear graduation. These tubes have consistent wall thickness with tight fit attached caps which ensures leak proof operation.

Catalog No.	Capacity (ml)
JW05631	0.5
JW05632	1.5
JW05633	2

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.