



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Centrifuge Tube Box

Product Image



Description

Material : Polypropylene

These Centrifuge tube boxes are molded in polypropylene . These boxes provide excellent substitute for conventional cardboard or foam racks. The transparent top gives a clear view of the tubes placed in the box which are easily identifiable because of the indexed grid. These boxes are autoclavable and can withstand temperatures ranging from -80°C to +121°C.

| Catalog No. | Capacity |
|--------------------|----------------------|
| 101439-A | For 15ml - 36 Places |
| 101439-B | For 50ml - 16 Places |

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.