



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Centrifuge Tube Box

## Product Image



# Description

---

**Material :** Polypropylene

These Centrifuge tube boxes are molded in polypropylene . These boxes provide excellent substitute for conventional cardboard or foam racks. The transparent top gives a clear view of the tubes placed in the box which are easily identifiable because of the indexed grid. These boxes are autoclavable and can withstand temperatures ranging from -80°C to +121°C.

Catalog No.	Capacity
101439-A	For 15ml - 36 Places
101439-B	For 50ml - 16 Places

# Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.