



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Centrifuge Tube Conical Bottom

## Product Image

---



## Description

---

**Material:** Polypropylene

Centrifuge tubes with conical bottom. Moulded in Polypropylene, these tubes are clear strong & safe. These autoclaveable tubes are provided with printed graduations on the outer surface. The built in seal ring in the cap makes these tubes absolutely leakproof. These tubes can with stand forces upto 7500xG12000G for 15 ml. tubes respectively. These tubes are available sterilized as well as non-sterilized.

<b>Catalog No.</b>	<b>Capacity (ml)</b>
JW05501	0.2
JW05502	0.5
JW05503	1.5
JW05504	15
JW05505	15 ml sterilized
JW05506	50
JW05507	50 ml sterilized

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.