



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Clement & Desormes's Apparatus

Product Image



Description

It consists of a glass reagent bottle of 10 litres capacity fitted with leak proof brass cap. The cap is fitted with one bent tube, one stopcock for connection to a pump and another stopcock for opening to air. It is provided with a football pump for filling air and one manometer on stand with scale. Without mercury.

Clement & Desorme's Apparatus

Catalog No. 100422

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.