



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Clinical Centrifuge, Economy

Product Image



Description

For 15ml. test tubes with fixed angle rotor head. It has 5 step speed regulator. Supply 230 volts, single phase 50Hz AC. Speed range from 0 to 3800 RPM.

Model	Catalog No.
Clinical Centrifuge, Economy, with fixed angle rotor head with 4 buckets	JA10140
Clinical Centrifuge, Economy, with fixed angle rotor head with 6 buckets	JA10141
Clinical Centrifuge, Economy, with fixed angle rotor head with 8 buckets	JA10142
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 4 buckets	JA10143
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 6 buckets	JA10144
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 8 buckets	JA10145

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.