



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Clinical Centrifuge, Economy

## Product Image



## Description

---

For 15ml. test tubes with fixed angle rotor head. It has 5 step speed regulator. Supply 230 volts, single phase 50Hz AC. Speed range from 0 to 3800 RPM.

Model	Catalog No.
Clinical Centrifuge, Economy, with fixed angle rotor head with 4 buckets	JA10140
Clinical Centrifuge, Economy, with fixed angle rotor head with 6 buckets	JA10141
Clinical Centrifuge, Economy, with fixed angle rotor head with 8 buckets	JA10142
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 4 buckets	JA10143
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 6 buckets	JA10144
Clinical Centrifuge, Economy, with fixed angle ABS rotor head with 8 buckets	JA10145

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.