



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Coagulation Analyzer

Product Image



Description

Catalog No. 100988

Coagulation Analyzer - Model 1232

TECHNICAL SPECIFICATIONS

Testing channels	4
Cuvette incubation positions	24
Reagent incubation positions	6
Wavelength	LED 470nm
Sample type	Plasma
Memory	10,000 results storage
Sample volume	20-40 μ L
Reagent volume	20-40 μ L
Electronic pipette	20-200 μ L Clotting factors are able to expressed in the form of time(Seconds), ratio, g/L Support LIS
Incubation temperature	37.0 \pm 1.0 $^{\circ}$ C
Relative humidity	\leq 85%
Deviation between channels	\leq 5%
Power supply	AC 100-240V 50/60Hz
Size	375mm x 335mm x 175mm
Net weight	5.2 Kg.

SALIENT FEATURES

- ✓ 5-inch touch screen and keypad 4 independent channels Parameters: PT, APTT, TT, FIB, etc.
- ✓ Principle: Clotting, optical method
- ✓ Micro volumes
- ✓ Disposable cuvette no carry-over
- ✓ Precise electronic pipette Languages: English, French, Spanish, Portuguese, etc.
- ✓ Automatic start at reagent addition
- ✓ Thermal printer, auto & manual print

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.