



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Coplin jar (Euro Design)

Product Image



Description

Material : Polypropylene

These Coplin Jar made of polypropylene can hold 10 microscope slides of the size 75mm x 25mm. Slides are held upright due to the grooved interior. The built in seal ring in the flat cap makes this Coplin Jar a completely leak proof one.

Catalog No.	Particulars
101440	Coplin jar (Euro Design)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.