



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cryo Box (PC)

Product Image



Description

Material: Polycarbonate

Cryo boxes are moulded in polycarbonate making it extra strong and autoclavable .These can accommodate cryp vials of up to 4.5 ml capacity and are designed to be used at temperature ranging from -190 to +120 ° C .The transparent cover of the bos gives a good view of the contents inside and is numbered for better inventory control.

Catalog No.	Particulars
JW07401	50 places for 1ml
JW07402	81 places for 1ml
JW07403	81 places for 4.5ml
JW07404	100 places for 1ml

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.