



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Cryo Box (PC)

## Product Image



## Description

---

**Material:** Polycarbonate

Cryo boxes are moulded in polycarbonate making it extra strong and autoclavable .These can accommodate cryp vials of up to 4.5 ml capacity and are designed to be used at temperature ranging from -190 to +120 ° C .The transparent cover of the bos gives a good view of the contents inside and is numbered for better inventory control.

<b>Catalog No.</b>	<b>Particulars</b>
JW07401	50 places for 1ml
JW07402	81 places for 1ml
JW07403	81 places for 4.5ml
JW07404	100 places for 1ml

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.