



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cryo Box (PP)

Product Image



Description

Material: Polypropylene

Moulded in polypropylene the top of the box is attached to the bottom through well designed built in hinges. The bottom has built in spaces to hold vials up to 2 ml capacity. Available in two different sizes. These autoclavable boxes bear alphanumeric index for easy identification of different vials placed inside.

Catalog No.

JW07501

JW07502

Particulars

81 places for 1ml

100 places for 1ml

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.