



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cryo Rack

Product Image



Description

Material: Polycarbonate

Cryo rack is moulded in polypropylene making it tough durable and autoclavable. This handy rack can hold 50 vials at a time. The unique feature of this rack is that vial can be unscrewed without actually picking them up as the vial will lock in every well and hence will not twist. This pack is provided with four anti skid rubbers and each well is identified with an alpha numeric index.

Catalog No. Particulars

JW07301 50 places for 1, 1.8 and 4.5ml Cryo Vials

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.