



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Cryo Rack

## Product Image



## Description

---

**Material:** Polycarbonate

Cyro rack is moulded in polypropylene making it tough durable and autoclavable. This handy rack can hold 50 vials at a time. The unique feature of this rack is that vial can be unscrewed without actually picking them up as the vial will lock in every well and hence will not twist. This pack is provided with four anti skid rubbers and each well is identified with an alpha numeric index.

**Catalog No. Particulars**

JW07301      50 places for 1, 1.8 and 4.5ml Cryo Vials

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.