



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cryo Vial

Product Image



Description

Material: Polypropylene

Cryo vial are normally used for dtroage of biologicalmaterial,human and animal cell.These sterilized tubes are provided with white prominent graduation and a long marking area for sample identification ,The cap of these vial is provided with a seal ring inside to make it completely leakproof.The longitudinal grovves in the base give a non twisting features to these vials when placed in a cryo racks.

Catalog No.	Capacity (ml)
JW07101	1
JW07102	1.8
JW07103	4.5

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.