



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Cyto Centrifuge, Brushless, 4000 R.P.M.

Product Image



Description

Model : 8

General Purpose Bench Top Laboratory Centrifuge Machine without Carbon:

- Heavy duty C.R.C. construction finished with powder coating.
- Fitted with brushless automatic balanced motor.
- Microprocessor based, 4 lines, 20 characters. LCD Display of set time, run time, set and run rpm, rcf. Equipped with 0-99 minute timer.
- Programmable speed regulator to control speed from 500 to 4000 R.P.M.
- 50 programmes can be stored in memory.
- Automatic rotor identification on display.
- Additional electronic lid lock which does not allow to open lid when machine is running.
- Dynamic brake (automatic) and imbalance detector are provided. Selection of 3 acceleration + 4 deceleration profile.

Technical Parameter :

Max. Speed	4000 r.p.m
Max. RCF	3600xg
Motor	Brushless Motor
Noise	585
Power	AC 220V 50 Hz 15A

Dimension:	(L x W x H) 490 x 380 x 270 mm
Weight:	29 kg

Cyto Centrifuge, Brushless, 4000 R.P.M

Catalog No. 100766

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.