



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Dental Glass Bead Sterilizer

## Product Image



## Description

Glass Bead Sterilizer is a very effective tool for sterilizing a wide range of Medical, Dental and specially endodontic instruments. It is portable in size and very convenient to use with automatic control to maintain desired temperature.

### Technical Specifications:

- Operating Voltage : 220V
- Size (H x W x D) : 14.5cm x 12.5cm x 13.5cm
- Well Size ( depth x dia.) : 6cm x 5cm

- Weight : 1.5 kg
- Quantity of glass beads: 250g/pc
- Start up time : 10-12 Min.
- Sterilization time : 10-20 seconds
- 

Dental Glass Bead Sterilizer

**Catalog No. 101003**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.