



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Digital Dissolved Oxygen Analyzer Cum Temperature Meter

## Product Image

---



## Description

---

Catalog No : 100118

Digital Dissolved Oxygen Analyzer Cum Temperature Meter - Model 27

## TECHNICAL SPECIFICATIONS

### DIGITAL D.O METER

✓ Range	0 to 20.0 ppm
✓ Temperature	0 to 100 (in Model <b>M-27</b> )
✓ Resolution	0.1 ppm
✓ Accuracy	0.1 ppm $\pm$ 1 count
✓ Temp. Compensation	0° C to 50° C Manual
✓ Compensation	$\pm$ 3 %
✓ D.O. Sensor	Amperometric Gold/Silver membrane type.

### GENERAL

✓ Display	3 ½ digit LED display
✓ Power	230V A.C 50 HZ
✓ Weight	2 Kg

### STANDARD ACCESSORIES

✓ D O Electrode	1 No.
✓ Instruction Manual	1 No.
✓ D.O Membranes	1 Set.
✓ Temperature Probe	1 No. (in Model <b>M-27</b> )

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.