



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Digital Leak Test Apparatus

## Product Image

---



## Description

---

Catalog No : 100128

Digital Leak Test Apparatus - Model 75

Precise Control of Vacuum

## TECHNICAL SPECIFICATIONS

✓ Type	GMP Model, In-built Test Unit
✓ Time Control & Range	Digital, Programmable, 00.01 sec to 99min59sec
✓ Test Vacuum Range	Adjustable up-to 650 mmHg
✓ Vacuum Display	2.5" SS Vacuum Gauge
✓ Vacuum Control & Release	By Single Front Panel Knob
✓ Audio Alarm at completion	Yes
✓ Vacuum Pump	Oil free, MV 27" +- 5%, 1/4HP Motor
✓ Test Unit Size	200mm, 250mm & 300mm (any one) (Specify in Order)

## ACCESSORIES

- ✓ Programmable Time With Digital Display
- ✓ Auto Stop After Completion Of Set Time With Alarm
- ✓ Precise Control of Vacuum
- ✓ Heavy Duty Oil-free Vacuum Pump
- ✓ 1mm Thick SS-304 cabinet

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.