



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Double Rotary Tablet Press

Product Image



Description

A double rotating tableting machine is built to function effectively, compressing tablets while ensuring optimum cleanliness and accessibility. As a result, the Rotary compression machine's upkeep is easy.

It can generate around 200000 pills each hour. Additionally, the machines can deliver tremendous compression pressure. To create definite and simple tablets from powdered

material this double rotating tableting machine is used. It is a machine that is extensively used in a variety of industries such as pharmaceuticals, herbal medicine manufacturing, ceramics, and other businesses linked to this industry.

High Quality product for use in Pharmaceutical Industry and Laboratory. Following is made in compliance with International Quality Standards and are manufactured using high-grade components and materials.

-

Double Rotary Tablet Press, 27 Stations	Catalog No. 101278 (A)
Double Rotary Tablet Press, 29 Stations	Catalog No. 101278 (B)
Double Rotary Tablet Press, 35 Stations	Catalog No. 101278 (C)
Double Rotary Tablet Press, 45 Stations	Catalog No. 101278 (D)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.