



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Draining Rack, Metal Wire

## Product Image

---



## Description

---

For general laboratory use. It is made of thick metal wire for durability and ruggedness. It is moreover covered with plastic layer to avoid any rusting. With an arrangement for hanging on the wall. Available in size & configuration of your choice. Custom made sizes available on demand. Modification possible as per requirements.

Plastic Type Draining Rack available in Plastic Lab Products Section.

<b>Catalog No.</b>	<b>Particulars</b>
JD20420	Draining Rack, Metal Wire

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.