



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Electrical C Arm OT Table

Product Image



Description

Model 501

- | Four sectional radiolucent top.
- | Eccentrically positioned top for maximum access for movement of C-ARM image intensifier.
- | Base and cover made of heavy & hygienic stainless steel sheet.
- | Mobile on smooth running castors and floor locking system for stability.
- | Leg section and head section are interchangeable.
- | Most useful for general surgery and certain special surgical procedure.
- | Smooth and accurate positioning by remote control for up – down, trendleburg/reverse trendleburg lateral tilt and head side raising lowering.

Technical Specification

- | Overall length of table top: - 1950mm.
- | Overall width of table top: -500mm.
- | Up-down movement 770 to 1020mm.
- | Kidney Section Up: - 150mm.
- | Trendelenburg position: - 25o.
- | Trendelenburg & Reverse position: - 25o.
- | Lateral tilt – 20o (Both sides).
- | Flex/Reflex position – 30o/220o.
- | Leg section - 90 o down manual.

Remote Control System Offers

- | Up and down position.

Standard Accessories

- | Anesthetic screen.
- | Shoulder supports with pads.
- | Arm boards with pads.
- | Side support with pads.
- | Black rubberized mattress.

Electrical C Arm OT Table

Catalog No. 100527

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.