



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Electrical C Arm OT Table

## Product Image



## Description

Model 501

- | Four sectional radiolucent top.
- | Eccentrically positioned top for maximum access for movement of C-ARM image intensifier.
- | Base and cover made of heavy & hygienic stainless steel sheet.
- | Mobile on smooth running castors and floor locking system for stability.
- | Leg section and head section are interchangeable.
- | Most useful for general surgery and certain special surgical procedure.
- | Smooth and accurate positioning by remote control for up – down, trendleburg/reverse trendleburg lateral tilt and head side raising lowering.

#### Technical Specification

- | Overall length of table top: - 1950mm.
- | Overall width of table top: -500mm.
- | Up-down movement 770 to 1020mm.
- | Kidney Section Up: - 150mm.
- | Trendelenburg position: - 25o.
- | Trendelenburg & Reverse position: - 25o.
- | Lateral tilt – 20o (Both sides).
- | Flex/Reflex position – 30o/220o.
- | Leg section - 90 o down manual.

#### Remote Control System Offers

- | Up and down position.

#### Standard Accessories

- | Anesthetic screen.
- | Shoulder supports with pads.
- | Arm boards with pads.
- | Side support with pads.
- | Black rubberized mattress.

Electrical C Arm OT Table

**Catalog No. 100527**

## Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.