



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Electrophoresis Power Supply, Analog Fixed Current

Product Image



Description

Model : EP200

Specifications

- Analog Model with fixed current 250mA and also 5 types of fixed voltage means 50V, 100V, 150V, 200V and 250V.
- It is recommended for mini Horizontal and mini vertical Electrophoresis Apparatus.
- The voltage is selected thru a 5 step voltage regulator. Example if selector switch is at 1 position. The LED light will glow at 50V means fixed voltage of 50V is selected. Like wise other voltage can be selected.

Electrophoresis Power Supply, Analog Fixed Current

Catalog No. 100783

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.