



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Electrophoresis System, Horizontal, Mini

Product Image



Description

Model : GE1

Specifications:

- It is required in biotechnology, Molecular Biology applications in colleges, Universities, Research Institutes.
- It is meant for Nucleic Acid (DNA-RNA) characterization.
- These are also called Submarine Electro-phoresis Unit/Agrose Gel/Poly
- Acrylamide Gel Electrophoresis (Mini size)
- The unit is jointless, moulded of transparent Polycarbonate to avoid leakage.
- It is fitted with 2 Platinum Electrodes and supplied with transparent jointless safety lid which is also made of polycarbonate.
- It is completed with Buffer chamber, safety lid with cables.

Supply as above and completed with following accessories

▶ Comb 9 Wells - 10 Teeths	(Set of 2 Pcs.)
▶ Comb 11 Wells - 12 Teeths	(Set of 2 Pcs.)
▶ Comb 13 wells - 14 teeths	(Set of 2 Pcs.)
▶ Casting Tray, Size 10x7cm	(Set of 1 Pcs.)
▶ Casting Tray, Size 7x7cm	(Set of 1 Pcs.)
▶ Red and Black Cable	(1 Each)

-
Electrophoresis System, Horizontal, Mini

Catalog No. 100781

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.