



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Electrophoresis System, Vertical, Mini

Product Image



Description

Model : GE2

Specifications:

- It is required in biotechnology, Molecular Biology applications in colleges, Universities, Research Institutes.
- It is meant for protein characterization.
- The buffer chamber is jointless, moulded of transparent Polycarbonate to avoid leakage. The unit is based on Acrylic sheet fitted with 3 leveling screws. Complete with Acrylic cover.

Technical Parameter

▶ Dimensions (Outer)(LxWxH)	15x18x15 cm.
▶ Gel Size (cm)	8x7 (1Gel)
▶ Buffer Required (ml)	100 + 150
▶ Platinum Electrodes	2 Nos.

Supply complete with standard accessories

• Wires black and red - 1 each
• Glass plate-1 and 1 notched rectangular Total 2nos.
• One set spacers of 1.0mm,1.5mm - 4 pcs. each
• 4 sets of Nylon moulded combs of 7 wells & 10 wells in 0.1mm and 1.5mm thickness
• Silicon Gasket - 1 No.

Electrophoresis System, Vertical, Mini

Catalog No. 100782

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.