



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Electrophoresis System, Vertical, Mini

## Product Image

---



## Description

---

Model : GE2

### Specifications:

- It is required in biotechnology, Molecular Biology applications in colleges, Universities, Research Institutes.
- It is meant for protein characterization.
- The buffer chamber is jointless, moulded of transparent Polycarbonate to avoid leakage. The unit is based on Acrylic sheet fitted with 3 leveling screws. Complete with Acrylic cover.

### Technical Parameter

▶ Dimensions (Outer)(LxWxH)	15x18x15 cm.
▶ Gel Size (cm)	8x7 (1Gel)
▶ Buffer Required (ml)	100 + 150
▶ Platinum Electrodes	2 Nos.

### Supply complete with standard accessories

• Wires black and red - 1 each
• Glass plate-1 and 1 notched rectangular Total 2nos.
• One set spacers of 1.0mm,1.5mm - 4 pcs. each
• 4 sets of Nylon moulded combs of 7 wells & 10 wells in 0.1mm and 1.5mm thickness
• Silicon Gasket - 1 No.

-  
Electrophoresis System, Vertical, Mini

**Catalog No. 100782**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.