



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Float Rack

Product Image



Description

Material: Polypropylene This float rack moulded in Polypropylene is quiet compact, autoclaveable and can hold 16 micro centrifuge tubes of 1.5ml for incubation in water baths. Easy to load, these are essential in microbiology & immunology labs. The surface of the rack has raised alpha numeric index for easy identification as has a central rod handle for convenient lifting.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.