



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Fluid Bed Dryer

## Product Image



## Description

It is a type of equipment that is widely used in the pharmaceutical industry to minimize the moisture content of medicinal powders and granules. It is intended for the rapid and even drying of powders, crystals, and medium-sized pellets.

High Quality product for use in Pharmaceutical Industry and Laboratory. Following is made in compliance with International Quality Standards and are manufactured using high-grade components and materials.

-

Fluid Bed Dryer, 20L	<b>Catalog No. 101243 (A)</b>
Fluid Bed Dryer, 100L	<b>Catalog No. 101243 (B)</b>
Fluid Bed Dryer, 220L	<b>Catalog No. 101243 (C)</b>
Fluid Bed Dryer, 440L	<b>Catalog No. 101243 (D)</b>
Fluid Bed Dryer, 780L	<b>Catalog No. 101243 (E)</b>
Fluid Bed Dryer, 960L	<b>Catalog No. 101243 (F)</b>

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.