



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Fluorescent Microscope, Research

Product Image



Description

- Body:** Fluorescent Microscope with 30° angle Binocular Head.
- Stage:** Regular stage 132x142mm.
- Focussing:** Co-axial coarse and fine focussing for easy manipulation based on four stage reduction gear system & travelling on ball bearing guide ways with virtually no wear, tear and free from backlash. Tension adjustment ring on one side. Illumination: 100 Watt Mercury Bulb with power source.
- Condenser:** Two lens Abbe condenser N.A. 1.25 with filter support, movable by rack & pinion arrangement.

- Eyepiece:** Wide Field 10X, 20mm
- Objective:** Head Objectives as specified below.
- Box:** Supplied in a thermocol packing.

Optional Accessories:

3W LED cool white light source in place of halogen light source.

Model	Catalog No.
Fluorescent Microscope with two objectives 40X (FL), 100X (FL).	JM004601
Fluorescent Microscope with four objectives 4X (FL), 10X (FL), 40X (FL), 100X (FL).	JM004602

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.