




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Gas Tap, Euro, 2 Way 90°, for Pharma & Lab


Product Image




Description

2 Way 90' Taps

- Taps are made from high quality Brass Forging Grade as per IS 8737 2007 & Machining grade as per IS : 319 (I) material to avoid porosity & leakage.
- Taps are Powder Coated with Full Glossy Epoxy Fusion Bond & are available in Grey Colour.
- Taps are tested in Pneumatic Test at 250 PSI, in Hydraulic Test at 30 PSI to avoid leakage.
- Life Cycle Test : Upto 5,00,000 Cycles.
- Water Static Pressure Max : 14 Bar



2 Way 90' Tap



Gas Tap, Euro, 2 Way 90°, for Pharma & Lab

Catalog No. 101182

Following tap also available for Water, Air, Nitrogen and Vacuum.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.