



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Gas Tap, Euro, 4 Way for Pharma & Lab

## Product Image

---




## Description

---

#### 4 Way Taps

- Taps are made from highly quality Brass Forging Grade as per IS 8737 2007 & Machining grade as per IS : 319 (I) material to avoid porosity & leakage.
- Taps are Powder Coated with Full Glossy Epoxy Fusion Bond & are available in Grey Colour.
- Taps are tested in Pneumatic Test at 250 PSI, in Hydraulic Test at 30 PSI to avoid leakage.
- Life Cycle Test : Upto 5,00,000 Cycles.
- Water Static Pressure Max : 14 Bar



Gas Tap, Euro, 4 Way for Pharma & Lab

**Catalog No. 101184**

Following tap also used for Water, Air, Nitrogen and Vacuum.

## Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.