



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Gynecology Delivery Bed

Product Image



Description

Model DT107

- | Stainless steel tubular Frame work.
- | Adjustable head section.
- | Movable leg section on 75mm castors, 2 with breaks.
- | 50 mm thick PU foam Mattress.
- | Accessories: - Lithotomic rods, side Push Handle, SS IV Rod and SS bowl & 50 mm PU foam Mattress.
- | Size : 1830L x 610W x 800H mm.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.