



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Haematocrit Centrifuge, 12000 R.P.M

## Product Image

---



## Description

---

Model : 15BL

### Specifications:

- Haematocrit Centrifuge Machine, Max. Speed 12000 r.p.m.
- Compact body which is made of plastic injection moulded.
- Fitted with maintenance free induction motor drive.
- Digital speed control range.
- Speed readout: Direct reading, digital tachometer.
- Safety lid-lock lid.
- Possibility of interrupting the run at any time by & quot, stop key.
- Timer range at least from 1 up to 15 min. in increments of 1 min. or in continuous.
- Possible quick spins.
- Diagnostic program and alarm, indicators (i.e. imbalance, over speed, lid open).
- Compatible with standard capillaries.
- Line voltage 220-240AC, 50-60Hz and fitted with a 13 amp plug.
- Microprocessor based Digital timer of range 10-15 min.
- Preset speed controls thru microprocessor from 2000 to 12000 r.p.m. in interval of 100 rpm.

### Supply completed with following accessories:

- Micro Haematocrit rotor to accommodate 24 capillaries.
- 75mm x 1mm or 12 capillaries.
- Rotor cover made of polycarbonate.
- Heparinised capillary tubes 75mm long x 1mm bore(Pack of 100 Nos.).
- Non-Heparinised capillary tubes 75mm long x 1mm bore (Pack of 100 Nos.).
- Reading device.
- Sealing wax- 1 tray

-

Haematocrit Centrifuge, 12000 R.P.M

**Catalog No. 100753**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.