



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Haematocrite Centrifuge, Digital, 13000 r.p.m

Product Image



Description

Compact table top Haematocrite Centrifuge with see thru acrylic lid, Fitted with variable speed controller, digital speed indicator, digital timer 0-15 minutes. Maximum rpm approx.13000 with load, Digital display shows Set time/ Run time /rpm, Microprocessor based Digital timer of range 0-15 min with LCD display, Preset speed controls thru microprocessor from 2000 to 13000 rpm

Supply completed with following :

Micro Haematocrite rotor to accommodate 24 capillaries of 75mmx1mm or 12 capillaries.
Heparinised capillary tubes 75mm long x 1 mm bore(Pack of 100 Nos.)
Non-Heparinised capillary tubes 75mm long x 1 mm bore(Pack of 100 Nos.)
Reading device.
Sealing wax- 1 tray.

Model	Catalog No.
Haematocrite Centrifuge, Digital, 13000 r.p.m	JA10210
Haematocrite Centrifuge, Digital, 13000 rpm with Brushless Motor	JA10211

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.