



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Hand Centrifuge, Heavy

## Product Image



## Description

---

For 15ml. test tubes with swing out head and polyamide buckets to house test tubes. It has clamping device for fixing on table. It gives 1500 to 2500 revolutions per minute.

<b>Model</b>	<b>Catalog No.</b>
Hand Centrifuge, Heavy, Bush type, 2 buckets	JA10120
Hand Centrifuge, Heavy, Bush type, 4 buckets	JA10121
Hand Centrifuge, Heavy, Ball Bearing, 2 buckets	JA10122
Hand Centrifuge, Heavy, Ball Bearing, 4 buckets	JA10123

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.