



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Hand Microtome, Bench Type

## Product Image



## Description

---

Hand Microtome is used to cut thin sections of Microscopic materials, such as specimen parts of plants and animals. It is considered as the basic form of a microtome permitting microtomic cuts by hand in a quick and accurate way. The microtome is suitable for practical courses in microscopy, for field work during excursions, for school education and for a serious amateur of biology.

### Specifications :

1. A well constructed hand operated bench microtome can be securely attached to the edge of any bench, table or work surface with the help of fixing clamp.
2. The section platform has a smooth wide flat top and durable surface to guide the sectioning razor & to protect the hand.
3. The feed screw has fine and long threads so that the embedded material may be advanced very gradually to cut thin sections easily.
4. The Razor is made of surgical quality steel with 75mm cutting edge.

### Particulars

### Catalog No.

Hand Microtome, Bench Type

JM00430

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.