



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Hydraulic C Arm OT Table

Product Image



Description

Model 502

-

- | Radio translucent top in four sections for use of C arm image intensifier.
- | Compatible for orthopedic attachments.
- | Heavy duty smooth gear mechanism to avail all operating positions.
- | Durable heavy duty leak proof hydraulic pump with four pillars.
- | Eccentric top and hydraulic chamber adjusted at one side of base for maximum use of CARM image intensifier.
- | Changeable oil valve for up & down movement of hydraulic pump.
- | Inter changeable head and leg section.
- | One radio translucent hand board and one hand rest with S.S. top for transfusion.

Technical Specification

- | Overall length of table top: - 1950mm.
- | Overall width of table top: -500mm.
- | Up-down movement 770 to 1020mm.
- | Kidney Section Up: - 150mm.
- | Trendelenburg position: - 25o.
- | Trendelenburg & Reverse position: - 25o.
- | Lateral tilt – 20o (Both sides).
- | Flex/Reflex position – 30o/220o.
- | Leg section - 90 o down manual.

Standard Accessories

- | Anesthetic screen.
- | Shoulder supports with pads.
- | Arm boards with pads.
- | Side support with pads.
- | Black rubberized mattress.

-

Hydraulic C Arm OT Table

Catalog No. 100528

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.