



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Hydraulic Delivery Bed

## Product Image

---



## Description

---

Model DT105

-

- | Mild Steel tubular frame work.
- | Height adjustable with foot operated Hydraulic pump.
- | Three section Top with High density PU form Mattress.
- | Mounted on 12 cm castors, 2 with brakes.
- | Provision of IV rod and Basin under top.
- | Side rail on the head side for holding hands.
- | Pre treated and epoxy powder coated finishing.

#### Technical Specification

- I) Height : - 500mm to 800mm ( $\pm$  25mm).
- ii) Total Length 800 to 1500 mm ( $\pm$  25mm).
- iii) Body section : 800mm x 300mm.
- iv) Leg section : 700mm x 300mm.

-  
Hydraulic Delivery Bed

**Catalog No. 100545**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.