




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Hydrometer, Density

Product Image



Description

For Analytical Laboratories. Shot & Wax poised, Adjusted For use at 15 °C. Accuracy: \pm 0.0002

Type A Hydrometer, Density, Calibration 0.001, Max Length 270 mm

Type B Hydrometer, Density, Calibration 0.0005, Max Length 335 mm

**Type A: Hydrometer,
Density, Calibration
0.001, Max Length
270 mm**

Catalog No.	Specific Gravity Range/ °Be Range
--------------------	--

JG26201	0.650-0.700 g/ml
---------	------------------

JG26202	0.700-0.750 g/ml
---------	------------------

JG26203	0.750-0.800 g/ml
---------	------------------

JG26204	0.800-0.850 g/ml
---------	------------------

JG26205	0.850-0.900 g/ml
---------	------------------

JG26206	0.900-0.950 g/ml
---------	------------------

JG26207	0.950-1.000 g/ml
---------	------------------

**Type B: Hydrometer,
Density, Calibration
0.0005, Max Length
335 mm**

Catalog No.	Specific Gravity Range/ °Be Range
--------------------	--

JG26301	0.650-0.700 g/ml
---------	------------------

JG26302	0.700-0.750 g/ml
---------	------------------

JG26303	0.750-0.800 g/ml
---------	------------------

JG26304	0.800-0.850 g/ml
---------	------------------

JG26305	0.850-0.900 g/ml
JG26306	0.900-0.950 g/ml
JG26307	0.950-1.000 g/ml

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.