



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Immuno Electrophoresis Apparatus

Product Image



Description

Comprising of walled migration cell made of white Perspex, outer chamber is provided with inlets and outlets for circulation of cold water, inner chamber is the buffer chamber fitted with electrodes with a clear transparent top lid fitted with a handle, a special Perspex carrier is supplied as micro immunological Accessory, which permits simultaneously study

of upto eight 3"x1" standard agar coated micro slides. Complete with plug, terminals and cables Without power pack. Size 30x30x8cm.

Particulars**Catalog No.**

Immuno Electrophoresis
Apparatus

JA20970

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.