



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan

## Product Image



## Description

---

Incubators are suitable for variety of uses like drying of slides, paraffin embedding, tissue culture work incubation of anti-body test, microbiological determination, incubation of hydroxy steroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of  $\pm 1^\circ \text{C}$ .

**Construction :** All external cases are made of mild steel, often with strong seamless welding & covered with a thick layer of stoved epoxy polyester paint, which provide a strong, highly screech & heat resistant surface for a long lasting attractive finish. Internal surfaces (Including the doors) are manufactured from stainless steel with great care in order to eliminate sharp edges. The double walled chamber is filled with high grade insulation to prevent heat loss. They are put inside beads, placed at the bottom & in both the side ribs for uniform temperature. The door closes on to an asbestos gasket & fitted with a glass window. The shelves slide neatly into the fixed runners. The fixed ventilators are provided on both side of the wall. The incubator is mounted on four rubber feet to prevent slipping & this protects the bench surface. The control panel houses a main on/off switch indicator lamp & temperature setting knob. Oven is controlled with capillary type thermostat. Incubator is provided with P.I.D. based controller for more precise control. It is also fitted with fan to ensure uniform heating in the chamber.

Particulars	Catalog No.
Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan, 300 x 300 x 300mm, 250 watts, 27 litre, 1 shelf	JA22180
Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan, 355 x 355 x 355mm, 300 watts, 45 litre, 2 shelves	JA22181
Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan, 405 x 405 x 405mm, 300 watts, 65 litre, 2 shelves	JA22182
Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan, 455 x 455 x 455mm, 300 watts, 95 litre, 2 shelves	JA22183
Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Fan, 605 x 605 X 605mm, 500 watts, 220 litre, 3 shelves	JA22184

**Note:** Custom-made incubators to your specific requirements like size, material, control mechanism, configuration, shelves or special needs available on demand.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change

them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.