



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic

Product Image



Description

Incubators are suitable for variety of uses like drying of slides, paraffin embedding, tissue culture work incubation of anti-body test, microbiological determination, incubation of hydroxy steroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of +/- 1° C.

Construction : All external cases are made of mild steel, often with strong seamless welding & covered with a thick layer of stoved epoxy polyester paint, which provide a strong, highly screech & heat resistant surface for a long lasting attractive finish. Internal surfaces (Including the doors) are manufactured from stainless steel with great care in order to eliminate sharp edges. The double walled chamber is filled with high grade insulation to prevent heat loss. They are put inside beads, placed at the bottom & in both the side ribs for uniform temperature. The door closes on to an asbestos gasket & fitted with a glass window. The shelves slide neatly into the fixed runners. The fixed ventilators are provided on both side of the wall. The thermometer is provided on the top of incubator to check the temperature. The incubator is mounted on four rubber feet to prevent slipping & this protects the bench surface. The control panel houses a main on/off switch indicator lamp & temperature setting knob. Oven is controlled with capillary type thermostat.

Particulars

Catalog No.

| | |
|---|---------|
| Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic, 300 x 300 x 300mm, 250 watts, 27 litre, 1 shelve | JA22130 |
| Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic, 355 x 355 x 355mm, 300 watts, 45 litre, 2 shelves | JA22131 |
| Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic, 405 x 405 x 405mm, 300 watts, 65 litre, 2 shelves | JA22132 |
| Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic, 455 x 455 x 455mm, 300 watts, 95 litre, 2 shelves | JA22133 |
| Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic, 605 x 605 X 605mm, 500 watts, 220 litre, 3 shelves | JA22134 |

Note: Custom-made incubators to your specific requirements like size, material, control mechanism, configuration, shelves or special needs available on demand.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.