



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan

Product Image



Description

Incubators are suitable for variety of uses like drying of slides, paraffin embedding, tissue culture work incubation of anti-body test, microbiological determination, incubation of hydroxy steroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of +/- 1° C.

Construction : External case is made of 304 grade stainless steel, with strong seamless welding & covered with a thick layer of stoved epoxy polyester paint, which provide a strong, highly screech & heat resistant surface for a long lasting attractive finish. Internal surfaces (Including the doors) are manufactured from stainless steel with great care in order to eliminate sharp edges. The double walled chamber is filled with high grade insulation to prevent heat loss. They are put inside beads, placed at the bottom & in both the side ribs for uniform temperature. The door closes on to an asbestos gasket & fitted with a glass window. The shelves slide neatly into the fixed runners. The fixed ventilators are provided on both side of the wall. The incubator is mounted on four rubber feet to prevent slipping & this protects the bench surface. The control panel houses a main on/off switch indicator lamp & temperature setting knob. Oven is controlled with capillary type thermostat. Incubator is provided with digital temperature controller cum indicator. It is also fitted with fan to ensure uniform heating in the chamber.

Particulars	Catalog No.
Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan, 300 x 300 x 300mm, 250 watts, 27 litre, 1 shelf	JA22170
Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan, 355 x 355 x 355mm, 300 watts, 45 litre, 2 shelves	JA22171
Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan, 405 x 405 x 405mm, 300 watts, 65 litre, 2 shelves	JA22172
Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan, 455 x 455 x 455mm, 300 watts, 95 litre, 2 shelves	JA22173
Incubator, Bacteriological, Memmert Type, Stainless Steel Inside and Outside, Digital Temperature Controller with Fan, 605 x 605 X 605mm, 500 watts, 220 litre, 3 shelves	JA22174

Note: Custom-made incubators to your specific requirements like size, material, control mechanism, configuration, shelves or special needs available on demand.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may

be made to design based on latest availability, process and design.