



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Incubator Thermometers

Product Image



Description

Pocket thermometer & Incubator Thermometers are specially designed by our creative designers and fabricated in compliance with industrial quality standards. These are available in perfect finish and provides perfect results. It is filled with mercury/ red toluene and are available with yellow robust plastic case with pen type clip.

Standard Available Ranges :

1- 00 to 500 C

2- 800 to 1200 F

Catalog No.

JD14925

Particulars

Incubator Thermometers

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.