



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Inverted Tissue Culture Microscope

Product Image



Description

- Body:** Binocular/Trinocular head
- Stage:** Triple layer co-axial mechanical stage for X-Y movement with fully graduated vernier scale for manipulation of specimen slide.
- Focussing:** Co-axial coarse and fine focussing for easy manipulation based on four stage reduction gear system & travelling on ball bearing guide ways with virtually no wear, tear and free from backlash. Tension adjustment ring on one side.
- Illumination:** 12V/50W incident lamp with variable light control arrangement.

Eyepiece: Wide Field 10X (Paired).

Objective: Long working distance objectives 5X, 10X, 20X, 40X SL.

Box: Supplied in wooden packing.

Model	Catalog No.
--------------	--------------------

Inverted Tissue Culture Microscope	JM005800
------------------------------------	----------

Note: Available with Trinocular head

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.