



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Laboratory Jack

## Product Image



## Description

An all purpose constructed from stainless steel, aluminum or mild steel. suitable for hot plates, burners, baths, glassware & many other items. It is a knob operated scissor action lifting mechanism for quick, easy & accurate height adjustment. Aluminum & mild steel jacks are finished with a thick layer of epoxy based polyester paint, which provides a strong highly scratch & chemical resistant surface.

Material	Height	Top Plate Size
----------	--------	----------------

Mild Steel	230 mm	15x15 cm
Mild Steel	330 mm	20x20 cm
Mild Steel	400 mm	25x25 cm
Stainless Steel 304 grade	230 mm	15x15 cm
Stainless Steel 304 grade	330 mm	20x20 cm
Stainless Steel 304 grade	400 mm	25x25 cm

Laboratory Jack

Catalog No. 100233

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.